

Ringbio

food safety solutions for dairy industry



heart of love

RINGBIO Rapid Tests ensure **your** milk safety





Ringbio® Rapid Tests for dairy milk analysis

100000	NBReader for Rapid Test Kit
100001	Beta Lactams Rapid Test
100002	Beta-Lactams + Tetracyclines, BT Combo Test
100003	Beta-Lactamase Rapid Test
100004	Aflatoxin M1 Rapid Test (0.05/0.1/0.2/0.3/0.4/0.5ppb)
100014	Chloramphenicol Rapid Test
100015	Beta/Tetra/Sulfa, 3-in-1 BTS TriTest
100016	Beta/Tetra/AflaM1, 3-in-1 TriTest BTM1 TriTest
100017	Beta/Tetra/Macrolides, 3-in-1 BTM TriTest
100018	Beta/Tetra/CAP/Strep, 4-in-1 BTCS QuaTest
100019	Beta/Tetra/Cefalexine, 3-in-1 BTCex TriTest
100020	Linco/Quino/Tylo, LQT 3-in-1 TriTest
100021	Beta-Lactams + Cefalexin, BCex Combo Test
100023	Macrolides Rapid Test, Tylosin/Tilmicosin 2-in-1 Test

NBReader® Rapid Test Strip Reader

Weight	3.0kg
Dimension	320mm × 195mm × 156mm
Wavelength	525 ± 30nm
Precision	C.V. ≤ 3%
Display	6.2Inch LCD
Printer	Built-in thermal printer
Memory	4GB built-in / flash drive

